
INSTRUBAC 20L

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
Cleaning protocol of Instrubac

CLEANING

Instrubac 20L is designed so as to be cleaned in an **automatic washer machine** at **max temperature 90°C** (Machines to comply with ISO15883).
Instrubac 20L can be manually cleaned using suitable detergent or disinfectant (see protocol established in your department).

- 1 Separate the basket from the tank.
- 2 Position the components separately in the washer machine.
Make sure the lid of the tank is open.

STERILISATION

Instrubac 20L can be **sterilised in a steam autoclave** to a **temperature 134° during 18 minutes** (max temperature : 137°C) (Autoclaves to comply with ISO17665).
The packaging materials and methods as well as their positioning in the autoclave must be taken into consideration to avoid deformation.

- 1 Place the basket into the tank. Keep the lid of the tank open.
- 2 Package up the whole thing with a sterilisation bag. Make sure the lid stays open.
- 3 Place the tank in the autoclave. Avoid all contact with other materials and products.

STORAGE

Untreated tanks must be stored in their original packaging.
Secondary packaging is recommended, for example, a cardboard box.
The tanks must be stored **in a dry place** (minimum 4°C) **without direct contact** with the ground.

REFERENCES

Simple box without tap.....Ref. 404 014
Double box without tap.....Ref. 404 019

Simple box with tap.....Ref. 404 021
Double box with tap.....Ref. 404 022

RECOMMENDATIONS :

- On first use, perform a manual sterilisation with adapted disinfectant product.
- A visual inspection is necessary to evaluate wear to the tank due to its use.
If it is damaged, do not use the tank and report it.