

Ref.: 80116
Código Sanitario: pendiente



Carga máxima 136 kg
Maximum load 136 kg



W4 es un andador práctico y ligero, de diseño minimalista, que te permitirá disfrutar de mayor autonomía en tus actividades diarias. Cuenta con manillar ajustable en altura para proporcionar un buen apoyo durante el paseo. Y cuando necesites descansar, su asiento acolchado de altura regulable te ofrece el mejor confort. Su plegado compacto y su ligereza hacen que sea muy fácil de almacenar y transportar. Y gracias a su gran cesta inferior y a la bolsa incorporada bajo el asiento, es muy práctico para llevar la compra y objetos personales sin necesidad de cargar con bolsas.

W4 is a practical and lightweight walker with a minimalist design that will allow you to enjoy greater autonomy in your daily activities. It features an adjustable handlebar height to provide good support during walking. When you need to rest, its padded and adjustable seat offers you the best comfort. Its compact folding and lightweight make it very easy to store and transport. Its large bottom basket and the built-in bag under the seat, it is very practical to carry shopping and personal items without the need to carry bags



Especificaciones Specifications

Medidas - Ancho, Largo, Alto Outer dimensions - Width, Long, Height	62,5 x 66,5 x 90 ~ 97 cm
Medidas plegado - Ancho, Largo, Alto Folding measures - Width, Long, Height	62,5 x 32 x 85 cm
Peso Weight	6,8 Kg
Max. Peso usuario Max. user weight	136 Kg
Homologación Standards	ISO 11199-2:2021

Medidas embalaje - Ancho, Largo, Alto Packaging measures - Width, Long, Height	1
Peso embalaje bruto Packaging gross weight	8 Kg
Medidas embalaje Master - Ancho, Largo, Alto Master Packaging measures - Width, Long, Height	57 x 57 x 16 cm
Peso embalaje Master bruto Master Packaging gross weight	19,25 kg
Unidades embalaje Master Units per Master package	1

Bolsa Desmontable Removable bag	✓
Medidas Bolsa Bag dimension	30 x 14 x 45 cm
Capacidad Bolsa Bag Capacity	19 L
Material del Tejido Fabric Material	Polyester
Tejido Hidrófugo Water Repellent Fabric	✓
Lavable a Máquina Machine Washable	✓
Cesta portaobjetos (desmontable) Storage basket (detachable)	✓

Acabado del Manillar Handlebar Finishing	Goma Rubber
Posiciones de manillar Handlebar positions	90 ~ 97 cm
Ruedas frontales Front wheels	6' (15,24 cm)
Bloqueo Ruedas frontales Front wheels lock	90° ~ 180° ~ 360°
Ruedas traseras Rear wheels	6' (15,24 cm)
Material del neumático Tyre material	ER ER
Freno de mano Hand Brake	✓
Reflectantes de seguridad Safety reflective	✓

Soporte para muleta Crutch holder	Accesorio Accessory
Posavasos Cup holder	Accesorio Accessory



Ref.: 80116 287